pNAPEP-9101

Substrat du tPA

(identique au Pefachrom tPA™ Pentapharm)

REF 61039101

RUO

Usage Recherche

CRYOPEP

83 rue Yves Montand
F-34080 MONTPELLIER
Tél.: +33(0)4 67 10 71 20
Fax : +33(0)4 67 10 71 21
contact@cryopep.com
www.cryopep.fr

COMPOSITION

Flacon de 25 mg de pNAPEP-9101, soit 0.0389 mmol

Excipient : (D)-mannitol (80 mg/flacon)

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Vial : 25 mg pNAPEP-9101, 0.0389 mmol

Bulking agent : (D)-mannitol (80 mg/vial)

The accurate quantity is indicated on analytical data sheet.

APPLICATION

This chromogenic substrate has been used for the determination of tissue plasminogen activator. It shows different sensitivity to native single-chain tPA (sc-tPA) and its fully active two-chain form (tc-tPA). This allows a differential determination of each molecular form.

REAGENT PREPARATION

Reconstitute the vial with 20 mL of sterile distilled water to obtain a concentration up to 1.945 mmol/L.

Let stand at room temperature (18-25°C) for 30 minutes.

Shake before use.

CHEMICAL STRUCTURE

Formula : CH₃SO₂-D-CHA-Gly-Arg-pNA, AcOH (Methanesulfonyl-D-cyclohexylalanin-glycyl-L-arginine-paranitroanilin acetate)

Molecular weight : 642.7 g/mol

PHYSICO-CHEMICAL DATA

Solubility : > 40 mmol/L (H₂O)

ε_{316 nm} : 1.27 . 10⁴ mole⁻¹ . L . cm⁻¹

STABILITY

Lyophilized powder :

Stable until the expiry date, if stored at 2-8°C. Avoid light exposure.

Keep in a dry place, hygroscopic product.

Aqueous solution (1 – 10 mmol/L) :

Stable for 3 months at 2-8°C.

PRECAUTIONS

- The product should be sealed in its original vial, following each use to ensure its stability.
- Minimize the evaporation of product during use.
- Contamination by microorganism may cause hydrolysis.
- Yellowing of the substrate indicates contamination, discard and use a new vial.

PRINCIPLE

The rate of pNA, released by enzymatic hydrolysis, is quantified by spectrophotometry at 405 nm and is proportional to the enzyme activity

ENZYMATIC DATA

The following kinetic constants were estimated at 37°C in Tris-imidazole buffer 50 mM pH 8.4, I 0.3.

Enzym	K _m (mmole/L)	V _{max} (nmole/µg tPA/min)
sc-t-PA one-chain	0.286	6.95
tc-t-PA two-chain	3 . 10 ⁻⁴	33.9

The tissue plasminogen activator is not affected by the following inhibitors used at the concentrations given.

Trasylol :	30 KIU/mL
Soybean trypsin inhibitor :	50 mg/mL
Antithrombin :	0.1 PEU/mL
Heparin :	3 UI/mL