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REFERENCES

Prasa D, Stuerzebecher J : Determination of activated factor IX in factor IX concentrates with a chromogenic substrate. Throm Res; 92:99-102

Stürzebecher J, Kopetzki E, Bode W, Hopfner KP : Dramatic enhancement of the catalytic activity of coagulation

Let stand at room temperature (18-25°C) for 30 minutes.

Shake before use.

CHEMICAL STRUCTURE

C23H36N8O7S, AcOH
Formula : CH3SO2-D-CHG-Gly-Arg-pNA, AcOH
(Methylsulfonyl-(D)-cyclohexylglycyl-glycyl-arginine-paranitroaniline monoacetate)
Molecular weight : 628,70 g/mol

PHYSICO-CHEMICAL DATA

Solubility : < 20 mmol/L (H₂O)
ε_{316 nm} : 1,27 . 10⁴ mole⁻¹ . L . cm⁻¹
ε_{405 nm} : 1 . 10⁴ mole⁻¹ . L . cm⁻¹

STABILITY


Lyophilized powder :
Stable until the expiry date, if stored at 2-8°C.
Avoid light exposure.
Keep in a dry place, hygroscopic product.
Shipment of product does not require cooling during the time of transportation.

Aqueous solution (1 - 4 mmol/L) :
Stable for 6 months at 2-8°C.

PRECAUTIONS

- The product should be sealed in its original vial, following each use to ensure its stability.
- Minimize the evaporation of product during use.
- Contamination by microorganism may cause hydrolysis.
- Yellowing of the substrate indicates contamination, discard and use a new vial.

PRINCIPLE

The rate of pNA, released by enzymatic hydrolysis, is quantified by spectrophotometry at 405 nm and is proportional to the enzyme activity

ENZYMATIC DATA

The following kinetic constants were estimated at 37°C in Tris buffer 50mM, pH 7.4, 100mM NaCl, 5mM CaCl₂, 40% (vol/vol) ethylene glycole.

Enzyme	K _m (mole/L)	k _{cat} (sec ⁻¹)
Factor IXa	1.3 . 10 ⁻³	4.4

REFERENCES

Prasa D, Stuerzebecher J : Determination of activated factor IX in factor IX concentrates with a chromogenic substrate. Throm Res; 92:99-102

Stürzebecher J, Kopetzki E, Bode W, Hopfner KP : Dramatic enhancement of the catalytic activity of coagulation

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COMPOSITION

Vial : 25 mg pNAPEP-9502, 0.0397 mmol
Bulking agent: (D)-mannitol (80 mg/vial)

The accurate quantity is indicated on analytical data sheet.

APPLICATIONS

This substrate could be used for the determination of activated factor IX. Its sensitivity is significantly enhanced in the presence of alcohols, especially ethylene glycol.

REAGENT PREPARATION

Reconstitute the vial with 20 mL of sterile distilled water to obtain a concentration around 2 mmol/L.