1<sup>ER</sup> / 3  
JUN  
2022

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MHEMO

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## CLINICAL VALIDATION OF THE T-TAS 01 PL ASSAY IN PATIENTS TAKING ANTIPLATELET THERAPY

J. Dahlen<sup>1</sup>, P. Gurbel<sup>2</sup>, A.H.B. Wu<sup>3</sup>, K. Bliden<sup>2</sup>, C-M. Ong<sup>3</sup>, U. Tantry<sup>2</sup>, K. Hosokawa<sup>4</sup>

<sup>1</sup>Hikari Dx, Inc., San Diego, United States, <sup>2</sup>Sinai Hospital, Baltimore, United States,

<sup>3</sup>University of California, San Francisco, United States, <sup>4</sup>Fujimori Kogyo Co., Ltd, Tokyo, Japan


Group	N	Mean	SD	Median	Range
Healthy Donors	142	381.5	55.5	390.9	142.5 – 467.7
AAS Monotherapy	57	218.4	114.4	205.7	2.7 – 410.9
Clopidogrel + AAS	18	46.2	47.3	21.7	3.6 – 159.8
Prasugrel + AAS	15	31.1	26.7	27.1	3.6 – 100.2
Ticagrelor + AAS	14	23.1	25.1	13.6	3.2 – 86.6

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# Comparaison T-TAS<sup>®</sup> 01 \_ PL Chip versus VerifyNow

Journal of Thrombosis and Haemostasis, 14: 850-859

DOI: 10.1111/jth.13256

## ORIGINAL ARTICLE

### Assessment of platelet-derived thrombogenicity with the total thrombus-formation analysis system in coronary artery disease patients receiving antiplatelet therapy

Y. ARIMA,\* K. KAIKITA,\* M. ISHII,\* M. ITO,\* D. SUETA,\* Y. OIMATSU,\* K. SAKAMOTO,\* K. TSUJITA,\* S. KOJIMA,\* K. NAKAGAWA,† S. HOKIMOTO\* and H. OGAWA\*

\*Department of Cardiovascular Medicine, Graduate School of Medical Sciences, Kumamoto University, Kumamoto, Japan; and †Division of Pharmacology and Therapeutics, Graduate School of Medical and Pharmaceutical Sciences, Kumamoto University, Kumamoto, Japan

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

# Exemple comparatif PL & AR chip

## A microchip flow-chamber assay screens congenital primary hemostasis disorders


Yuto Nakajima, MD,<sup>1</sup> Koji Yada, MD, PhD,<sup>1,2</sup> Kenichi Ogiwara, MD, PhD,<sup>1</sup> Shoko Furukawa, MD, PhD,<sup>1,3</sup> Naruto Shimonishi, MD,<sup>1</sup> Midori Shima, MD, PhD,<sup>1</sup> Keiji Nogami, MD, PhD<sup>1</sup>

**Methods:** Microchips coated with collagen (PL-chip) and collagen/thromboplastin (AR-chip) were utilized to evaluate platelet thrombus formation (PTF) at high shear and fibrin-rich PTF at low shear, respectively, in whole bloods from 22 VWD (16 type 2A, 4 type 2B, 2 type 3) and 4 PFDs (2 Bernard-Soulier syndrome, 2 Glanzmann's thrombasthenia) patients. The time-to-increase by 10 kPa ( $T_{10}$ ) was calculated from flow pressure curves. Also, whole blood-induced platelet aggregation was assessed by Multiplate<sup>®</sup>.

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


**Fig. 2** T-TAS discriminates hemostatic function of major bleeding patients from that of minor bleeding patients. Overlays of T-TAS waveforms of 10 patients requiring platelet (PLT) transfusion, upper panel. Blue and red lines indicate waveforms in the minor bleeding ( $n = 4$ ) and major bleeding ( $n = 6$ ) patients, respectively. The area under the curve (AUC) of the T-TAS, expressed as mean  $\pm$  standard deviation, lower panel. The significance of the difference between the major and minor bleeding groups was analyzed with the Student  $t$  test. **\*\*** $P < 0.01$

RESEARCH

Open Access

## A modified microchip-based flow chamber system for evaluating thrombogenicity in patients with thrombocytopenia

Bengo Atari<sup>1</sup>, Takashi Ito<sup>2\*</sup>, Tomoka Nagasato<sup>3</sup>, Tomoko Ohnishi<sup>3</sup>, Kazuya Hosokawa<sup>3</sup>, Tomotsugu Yasuda<sup>1</sup>, Ikuro Maruyama<sup>2</sup> and Yasuyuki Kakihana<sup>1</sup>

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# T-TAS<sup>®</sup> 01

Total Thrombus formation Analysis System

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Can the Total Thrombus-formation Analysis System (T-TAS) better predict coagulation disorders than conventional laboratory measurements in patients with polycystic ovary syndrome?

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**Endometriosis is associated with an increased whole-blood thrombogenicity detected by a novel automated microchip flow-chamber system (T-TAS<sup>®</sup>)**

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- ESC Août, 2022

- ✓ Cardiovascular / Antithrombotic therapy

Accès aux publications

<https://www.t-tas.info/pub/>

## ESC Congress 2022 Barcelona

ONSITE & ONLINE,  
26-29 AUGUST

THE MAGIC  
OF CARDIOLOGY  
#ESCCongress