pNAPEP-1703

Substrat de la plasmine
(identique au S-2403™ de Chromogenix)

REF 61011703

RUO

Usage Recherche

CRYOPEP

83 rue Yves Montand
F-34080 MONTPELLIER
Tél.: +33(0)4 67 10 71 20
Fax : +33(0)4 67 10 71 21
contact@cryopep.com
www.cryopep.fr

COMPOSITION

Flacon de 25 mg de pNAPEP-1703, soit 0.0446 mmol

Excipient : (D)-mannitol (80 mg/flacon)

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Vial : 25 mg pNAPEP-1703, 0.0446 mmol

Bulking agent : (D)-mannitol (80 mg/vial)

The accurate quantity is indicated on analytical data sheet.

APPLICATION

This chromogenic substrate could be used for the determination in plasma of :

- antiplasmin
- plasminogen
- activity plasmin-like

REAGENT PREPARATION

Reconstitute the vial with 20 mL of sterile distilled water to obtain a concentration up to 2.23 mmol/L.

Let stand at room temperature (18-25°C) for 30 minutes.

Shake before use.

CHEMICAL STRUCTURE

C₂₆H₃₂N₆O₆, HCl

Formula : pGlu-Phe-Lys-pNA, HCl
(L-Pyroglutamyl-L-Phénylalanyl-L-Lysine-paranitroaniline hydrochloride)

Molecular weight : 561.0 g/mol

PHYSICO-CHEMICAL DATA

Solubility : > 40 mmol/L (H₂O)

$\epsilon_{316\text{ nm}}$: $1,27 \cdot 10^4 \text{ mole}^{-1} \cdot \text{L} \cdot \text{cm}^{-1}$

STABILITY

Lyophilized powder :

Stable until the expiry date, if stored at 2-8°C.

Avoid light exposure.

Keep in a dry place, hygroscopic product.


Aqueous solution (3 - 8 mmol/L) :

Stable for 6 months at 2-8°C.

PRECAUTIONS

- The product should be sealed in its original vial, following each use to ensure its stability.
- Minimize the evaporation of product during use.
- Contamination by microorganism may cause hydrolysis.
- Yellowing of the substrate indicates contamination, discard and use a new vial.

PRINCIPLE

The rate of pNA, released by enzymatic hydrolysis, is quantified by spectrophotometry at 405 nm and is proportional to the enzyme activity

ENZYMATIC DATA

The following kinetic constants were estimated at 37°C in Tris buffer pH 8.4, I 0.15.

Enzyme

Human plasmin : $K_m = 2.9 \times 10^{-4} \text{ mol/L}$; $k_{cat} = 92 \text{ s}^{-1}$

Plasminogen SK : $K_m = 3.0 \times 10^{-4} \text{ mol/L}$; $k_{cat} = 43 \text{ s}^{-1}$

SELECTIVITY

This substrate is insensitive to plasma-kallikrein, thrombin and factor Xa.