pNAPEP-8703

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Vial : 25 mg pNAPEP-8703, 34.92 µmol

Bulking agent : (D)-mannitol (80 mg/vial)

The accurate quantity is indicated on analytical data sheet.

APPLICATIONS

The chromogenic substrate pNAPEP 8703 is highly sensitive for C1-esterase and is used for the determination of C1-esterase inhibitor (C1 INH) in plasma. C1 INH is a regulatory protein that acts as an inhibitor of various serine proteases in the complement system, the kallikrein-kinin system, the coagulation cascade and in fibrinolysis.

REAGENT PREPARATION

Reconstitute the vial with 17,5 mL of sterile distilled water to obtain a concentration around 2.0 mmol/L.

Let stand at room temperature (18-25°C) for 30 minutes.

Shake before use.

CHEMICAL STRUCTURE

Formula : CH₃CO-Lys(Cbo)-Gly-Arg-pNA, AcOH (Methylcarbonyl-lysyl(ε-benzyloxycarbonyl)-glycyl-arginine-paranitroaniline monoacetate)

Molecular weight : 715,80 g/mol

PHYSICO-CHEMICAL DATA

Solubility : >10 mmol/L (H₂O)

ε_{316 nm} : 1.27 . 10⁴ mole⁻¹ . L . cm⁻¹

ε_{405 nm} : 1 . 10⁴ mole⁻¹ . L . cm⁻¹

STABILITY

Lyophilized powder :

Stable until the expiry date, if stored at 2-8°C in the dark. Avoid light exposure.

Keep in a dry place, hygroscopic product.

Shipment of product does not require cooling during the time of transportation.

Aqueous solution (1 - 5 mmol/L) :

Stable for 6 months at 2-8°C.

PRECAUTIONS

- The product should be sealed in its original vial, following each use to ensure its stability.
- Minimize the evaporation of product during use.
- Contamination by microorganism may cause hydrolysis.
- Yellowing of the substrate indicates contamination, discard and use a new vial.

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The rate of pNA, released by enzymatic hydrolysis, is quantified by spectrophotometry at 405 nm and is proportional to the enzyme activity.

ENZYMATIC DATA

Suggested protocol :

Plasma 100 µL

C1-estérase 100 µL

⇒ Incubate for 5 mn at 37°C

Tampon 700 µL

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