❖ THROMBOINCODE

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^(#) Homozygotic value

^(#) Data obtained from literature and diverse meta-analyses

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
PANEL GENETIQUE

1

Contact activation (intrinsic) pathway

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Study conducted on the Spanish population including 249 cases and 248 controls
 * % of patients with VTD in whom at least one genetic risk factor for thrombosis is detected

Description pathologie

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Exemple

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FAQ

Bibliographie

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2

Description pathologie

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Exemple

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FAQ

Bibliographie

	VALUES	INDIVIDUAL	
Clinical Data	Gender	M=1 F=2	2
	Age	Years	35
	Smoker		0
	Type II Diabetes		0
	Family history	Yes=1	1
	Pregnancy	No=0	0
	Oral contraceptives		1
	BMI	Kg/m ²	25
Genetic Data	FV Leiden		0
	Prothrombin		0
	Serpin A10		0
	FXII	No=0	0
	FXIII	Hetero=1	0
	Serpin C1	Homo=2	1
	A1 Carrier		1
	FV Cambridge		0
	FV Hong Kong		0

Clinical-genetic risk of thrombosis^(*): 16.8%

* Clinical-genetic risk of thrombosis: 16.8 % means that 17 out of every 100 patients with the same clinical-genetic conditions will suffer VT within 10 years.

Algorithme : le score

2

Description pathologie

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Exemple

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FAQ

Bibliographie

	PATIENT A	PATIENT A
Gender Men=1 Women=0	0	0
Age	35	35
FV Leiden 0=no; 1-Hetero; 2-Homo; 3-Carrier	0	0
Prothrombin alone	0	0
Protein Z inh (SerpinA10)	0	
FXII	0	
FXIII	0	
AT Cambridge (SerpinC1)	1	
ABO: A1 carrier	1	
OC alone	1	1
Diabetes type 2	0	0
Smoking	0	0
Family history	1	1
Pregnacy	0	0
BMI	26	26
Calculated risk	16.8%	2.5%

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	PATIENT A	PATIENT A
Gender Men=1 Women=0	0	0
Age	35	35
FV Leiden 0=no; 1-Hetero; 2-Homo; 3-Carrier	0	0
Prothrombin alone	0	0
Protein Z inh (Serpin A10)	0	
FXII	0	
FXIII	0	
AT Cambridge (Serpin C1)	1	
ABO: A1 carrier	1	
OC alone	1	1
Diabetes type 2	0	0
Smoking	0	0
Family history	1	1
Pregnacy	0	0
BMI	26	26
Calculated risk	16.8%	2.5%

Description pathologie

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Exemple

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FAQ

Bibliographie

La synthèse

3

Description pathologie

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Exemple

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FAQ

Bibliographie

30001155; 2013-01-28
Results report, Page 2 of 2

Variant	Genotype	Interpretation
F12, 46C>T	C/C	Two normal alleles
ABO Blood group, rs8176719	delG/G	The patient carries one A1 allele at ABO locus
ABO Blood group, rs7853989	C/C	
ABO Blood group, rs8176743	G/G	
ABO Blood group, rs8176750	C/C	
SerpinA10, Arg67Stop	C/C	Two normal alleles (Arg67)
SerpinC1, Ala384Ser (Cambridge II)	G/G	Two normal alleles (Ala384)
F5, Arg506Gln (F5 Leiden)	G/G	Two normal alleles (Arg506)
F5, Arg306Thr (F5 Cambridge)	G/G	Two normal alleles
F5, Arg306Gly (F5 Hong Kong)	A/A	Two normal alleles
F13, Val34Leu	G/G	Two risk alleles (Val34)
Prothrombin, G20210A	G/G	Two normal alleles

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	PATIENT A
Gender Men=1 Women=0	0
Age	35
FV Leiden 0=no; 1=Hetero; 2=Homo; 3=Carrier	0
Prothrombin alone	0
Protein Z inh (SerpinA10)	0
FXII	0
FXIII	0
AT Cambridge (SerpinC1)	1
ABO: A1 carrier	1
OC alone	1
Diabetes type 2	0
Smoking	0
Family history	1
Pregnacy	0
BMI	26

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
2/4

THROMBOINCODE

Individualized Recommendation and Genetic Counselling Report

Sample code:	0000004	Sample type:	Saliva	Request:	22/03/2012
Lab Code.:	6784563	Requesting physician:	J. Smith	Delivery:	28/03/2012

Reason for consultation:

Patient with a history of venous Thromboembolism (VTE)

Person with a family history (especially first degree relatives) of VTE and who faces a situation of risk of venous thrombosis:

- medication: oral contraceptives

Comments: DVT episode after placing a portacath during chemotherapy for a multifocal infiltrating adenocarcinoma of the large intestine (no thrombophilia test performed).

Without other factors that counterindicate the use of oral contraception, we wish to know the patient's genetic risk of thrombosis to assess the use of oral contraceptives

Clinical Data:

Patient Data:

Current age: 25 years Gender: Female Weight: 59 kg Height: 172 cm

Blood pressure: 110 / 70 mmHg

Family history of venous thromboembolic disease (VTD):

Family history of VTD

Family history of artery thrombosis:

No family history of artery thrombosis

Personal cardiovascular risk factors:

No known history

Exemple

3/4

THROMBOINCODE

Individualized Recommendation and Genetic Counselling Report

Sample code:	0000004	Sample type:	Saliva	Request:	22/03/2012
Lab Code.:	0784563	Requesting physician:	J. Smith	Delivery:	28/03/2012

Dear Doctor:

We have used the Thrombo inCode platform to analyze the above mentioned sample and have tested for polymorphisms in variants associated with the risk of developing thrombosis.

Results of the genetic variant test:

Gene	Variant	Risk allele	Result	
			Genotype	Interpretation
ABO	rs8176719	na,**	n.a.	Carrier of an A1 allele at the ABO locus
ABO	rs7853989	na,**		
ABO	rs8176743	na,**		
ABO	rs8176750	na,**		
F12	46C>T	T	C/T	Carrier of heterozygous risk allele 46T
Serpin A10	Arg67Stop	T	C/C	No risk alleles have been identified
Serpin C1	Ala384Ser	T	G/G	No risk alleles have been identified
F5	Arg506Gln (F5 Leiden)	A	G/A	Heterozygous carrier of Leiden mutation (Gln506)
F5	Arg306Thr (F5 Cambridge)	C	G/G	No risk alleles have been identified
F5	Arg306Gly (F5 Hong Kong)	G	G/G	No risk alleles have been identified
F13	Val34Leu	T	T/T	No risk alleles have been identified
Protrombin	G20210A	A	G/G	No risk alleles have been identified

* Protecting allele

** na (not applicable). These alleles are assessed according to the ABO genotype.

Phenotype Interpretation:

Thrombo inCode estimates the overall risk of having a thrombotic event based on a validated algorithm that includes genetic and clinical data. According to this algorithm the risk for your patient of developing a thromboembolic event is 5,8%.

Heterozygous carriers of the Arg506Gln mutation (FV Leiden) in the factor V gene (F5) show an increase of a phenotype known as "Resistance to activated protein C, RPCA", and their risk of developing a venous thrombosis is 4-5 times greater than non-carriers.

Heterozygous carriers of the Arg506Gln mutation (FV Leiden) in the factor V gene (F5) and whose ABO group is non-OO (specially group A1) have 2-5 times more risk of developing a venous thrombosis than group OO.

General Recommendations / Observations:

According to the clinical data provided and the results of the genetic profile, the following should be considered for your patient:

- Anticoagulant prophylaxis in situations of thrombo-embolic risk.
- Avoid, if possible, the use of oral contraceptives and hormone replacement therapy.
- Avoid, if possible, the development of other risk factors.
- Family study (first-degree kinship: parents, children, siblings)

Description pathologie

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Exemple

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FAQ

Bibliographie

Exemple

Description
pathologie

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Exemple

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FAQ

Bibliographie

THROMBOINCODE

Individualized Recommendation and
Genetic Counselling Report

4/4

Sample code:	0000004	Sample type:	Saliva	Request:	22/03/2012
Lab Code.:	6784563	Requesting physician:	J. Smith	Delivery:	28/03/2012

Remark:

A person with the same characteristics but without a genotype associated with thrombosis would have a 0,2% risk for developing thrombosis. Thus, this patient currently has 29 times more risk of developing thrombosis than a person without the genetic background.

The use of oral contraception increases the risk for developing thromboembolic events. In this patient, oral contraception would increase the risk from 5,8 to 15,9%. She would have 34 times more predisposition than a person with her characteristics, but without the genotype associated to thrombosis.

The above-mentioned risks can significantly increase if the patient has other risk factors for thromboembolic events (smoking, hypertension, obesity, or diabetes mellitus). Therefore, we recommend controlling these other risk factors.

With regard to the reason for consultation, our recommendation is to avoid the use of oral contraception.

It is very likely that the thrombotic event that the patient's mother had is associated to a risk genotype.

Barcelona, March 28th 2012

Dr. José Manuel Soria
Head of Genetic Diagnosis Team

Literature:

1. Tirado I, Soria JM, Mateo J, Oliver A, Souto JC, Santamaria A, Felices R, Borrell M, Fontcuberta J. Association after linkage analysis indicates that homozygosity for the 4600 T polymorphism in the F12 gene is a genetic risk factor for venous thrombosis. *Thromb Haemost* 2004;91:896-904.
2. Corral J, González-Conejero R, Soria JM, González-Porras JR, Pérez-Ceballos E, Lecumberri R, Roldán V, Souto JC, Miñano A, Hernández-Espinosa D, Alberca I, Fontcuberta J, Vicente V. A Nonsense Polymorphism in the Protein Z-dependent Protease Inhibitor Increases the Risk for Venous Thrombosis. *Blood* 2008;108:177-183.
3. Corral J, Hernández-Espinosa D, Soria JM, González-Conejero R, Ordoñez A, González-Porras JR, Pérez-Ceballos E, Lecumberri R, Sanchez I, Roldán V, Mateo J, Miñano A, González M, Alberca I, Fontcuberta J, Vicente V. Antithrombin Cambridge II (Ala384Ser): an underestimated but strong and frequent genetic risk factor for venous thrombosis. *Blood* 2007;109:4258-4263.
4. Bezemer ID, Bare LA, Doggen CJM, et al. Gene variants associated with deep vein thrombosis. *JAMA* 2008;299:1306-1314.
5. Lindhoff-Last E, Luxembourg B. Evidence-based indications for thrombophilia screening. *Vasa* 2008;37:19-30.
6. Wells PS, Anderson JL, Scarvelis DK, Doucette SP, Gagnon F. Factor XIII Val 34Leu variant is protective against venous thromboembolism: A HuGE review and meta-analysis. *Am J Epidemiol* 2008;164:101-109.
7. Wu O, Bayoumi N, Vickers MA, Clark P. ABO(H) blood groups and the risk of vascular disease: a systematic review and meta-analysis. *J Thromb Haemost* 2008;6:62-69.
8. Chan WP, Lee CK, Kwong YL, Lam CK, Liang R. A novel mutation of Arg306 of factor V gene in Hong Kong Chinese. *Blood* 1998;91:1135-1139.
9. Williamson D, Brown K, Ludington R, Baglin T. Factor V Cambridge: a new mutation (Arg306→Thr) associated with resistance to activated protein C. *Blood* 1998;91:1140-1144.
10. Cochery-Nouvellon E, Mercier E, Lissalde-Lavigne G, Daurès JP, Quéré I, Dauzat M, Marès P, Gris JC. Homozygosity for the C46T polymorphism of the F12 gene is a risk factor for venous thrombosis during the first pregnancy. *J Thromb Haemost* 2007;5:700-707.
11. Reuner KH, Jenetzky E, Aleu A, Litfin F, Mellado P, Kloss M, Jüttler E, Grau AJ, Rickmann H, Patschke H, Lichy C. Factor XII C46T gene polymorphism and the risk of cerebral venous thrombosis. *Neurology* 2008;70:129-132.

The genetic study analyzes 12 variants in 7 genes associated with thrombotic risk. This service uses a microarray

Exemple

Description pathologie

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Exemple

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FAQ

Bibliographie

THROMBOINCODE PERSONALIZING RISK ASSESSMENT IN THROMBOSIS
APPLICATOR ASSESSM

Requesting Physician*:

Name*: Dr. (Surname) _____ (Name)	Specialty*:
Professional ID number*:	Contact phone number*:
e-mail address:	Country*:
Institution*:	

Patient identification*:

Patient Code:	Extraction date*: _____ (dd/mm/yyyy)
Report language <input type="radio"/> English <input type="radio"/> Spanish	Sample type: <input type="radio"/> Saliva <input type="radio"/> Blood <input type="radio"/> DNA <input type="radio"/> Other (specify) _____
Your reference: _____	

1-Reason for the consultation (more than one option can be selected)*:

Patient with a history of venous Thromboembolism (VTE). **Section 2 should be answered.**

Patient with history of arterial thrombosis (myocardial infarction, stroke or peripheral arterial disease). **Section 3 should be answered.**

Patients with history of: repeated abortions / implantation failures preeclampsia

Person with a family history (especially first degree relatives) of VTE and who faces a situation with risk of venous thrombosis:
 oral contraceptives other medication pregnancy surgical intervention other _____

Family of carriers of genetic variants that predispose to thrombosis _____

Other reasons. Please, specify _____

Comments: _____

2-Personal history of VTE*

No history

Deep venous Thrombosis (DVT)

Pulmonary embolism

Deep venous Thrombosis (DVT)

Single episode DVT recurrence

Location of the first episode:

Lower limb	Venous cerebral sinus	Retina
Upper limb	Visceral (Mesenteric, Splenoportal, Renal)	Recurrent superficial thrombophlebitis (STP)
Other		

Pulmonary embolism.

Isolated Secondary to DVT

Data of the first DVT episode

Spontaneous Secondary to:

Surgery	Traumatism
Immobilization	Oral contraceptives
Pregnancy /postpartum	Hormonal Replacement Therapy
Active Cancer	Other

Age at first episode _____ years

Number of different VTE episodes (PE+DVT+STP) _____

Exemple

Description pathologie

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Exemple

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FAQ

Bibliographie

THROMBOINCODE PERSONALIZING RISK-ASSESSMENT IN THROMBOSIS

2/2 APPLICATION ASSESSM

3.-Personal history of arterial thrombosis*

<input type="checkbox"/> No history			
<input type="checkbox"/> Myocardial infarction	Age at 1st		
<input type="checkbox"/> Coronary angina	Age at 1st		
<input type="checkbox"/> Ischemic Ictus (as T.I.A.)	Age at 1st		<input type="checkbox"/> Cardioembolic
<input type="checkbox"/> Ischemic Ictus (cerebral infarction)	Age at 1st		<input type="checkbox"/> Cardioembolic
<input type="checkbox"/> Peripheral arterial disease	Age at 1st		<input type="checkbox"/> Cardioembolic
<input type="checkbox"/> Bypass or arterial stent Thrombosis			
<input type="checkbox"/> Absence of arterial risk factors (hyperliprotein, dyslipidemia, obesity, smoking, diabetes), arteriosclerosis or other obvious cause of vascular ischemia			

4.-Patient's characteristics and personal risk factors*

Age years Gender: Female Male
 Weight Kg lbs. Height cm ft. in.

Risk factors

<input type="checkbox"/> Family History of Thromboembolic Disease:	<input type="checkbox"/> Venous Thrombosis	<input type="checkbox"/> Arterial Thrombosis
<input type="checkbox"/> Hypertension	<input type="checkbox"/> Smoking habits	<input type="checkbox"/> Diabetes
<input type="checkbox"/> Dyslipemia	<input type="checkbox"/> Active cancer	<input type="checkbox"/> Chemotherapy
<input type="checkbox"/> Neurologic disease with limbs paralists	<input type="checkbox"/> Venous catheter or pacemaker carrier	<input type="checkbox"/> Inflammatory bowel disease
<input type="checkbox"/> Other autoimmune diseases: <input type="text"/>		
<input type="checkbox"/> Prothrombotic drugs:	<input type="checkbox"/> Oral contraceptives	
	<input type="checkbox"/> Hormone replacement therapy	
	<input type="checkbox"/> Therapy with selective estrogen receptor modulator, SERM	
<input type="checkbox"/> Other: <input type="text"/>		
<input type="checkbox"/> Pregnancy or puerperium		

5.-In women:

Number of pregnancies:

Spontaneous abortion / Implantation failures: 1st quarter 2nd quarter 3rd quarter

Obstetric complications Eclampsia Preeclampsia

Intrauterine growth retardation

I declare that the patient has signed an Informed Consent for the genetic test.

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
Bibliographie

Bibliographie

- 1 Ariëns RA et al. The factor XIII V34L polymorphism accelerates thrombin activation of factor XIII and affects cross-linked fibrin structure. *Blood*. 2000; 96(3):988-95.
- 2 Bates SM et al. VTE, thrombophilia, antithrombotic therapy and pregnancy. *Antithrombotic Therapy and the Prevention of Thrombosis, 9th ed: American College of Chest Physicians Evidence-based Guidelines*. *Chest*. 2012 Feb; 141(2)(Suppl):e691S-e736S
- 3 Corral J et al. A nonsense polymorphism in the protein Z-dependent protease inhibitor increases the risk for venous thrombosis. *Blood*. 2006; 108(1):177-183.
- 4 Corral J et al. Antithrombin Cambridge II (A384S): an underestimated genetic risk factor for venous thrombosis. *Blood*. 2007; 109(10):4258-4263.401.
- 5 Cochery-Nouvellon E et al. Homozygosity for the C46T polymorphism of the F12 gene is a risk factor for venous thrombosis during the first pregnancy. *J Thromb Haemost*. 2007; 5(4):700-7.

Description
pathologie

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Exemple

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
FAQ

Bibliographie

Bibliographie

- 6 Larsen TB, Johnsen SP, Gislum M, Møller CAI, Larsen H, Sørensen HT. ABO blood groups and risk of venous thromboembolism during pregnancy and the puerperium. A population-based, nested case-control study. *J Thromb Haemost* 2005;3:300-304.
- 7 Norstrem E et al. Functional characterization of recombinant FV Hong Kong and FV Cambridge. *Blood* 2002; 100: 524-530.
- 8 Salas E et al. Predicting the risk of thromboembolism disease by analyzing multiple genetic variants: incorporating risk profiles into clinical practice. Eurothrombosis Lisboa, 7th-8th October 2011. Poster.
- 9 Schroeder V et al. New developments in the area of factor XIII. *J Thromb Haemost* 2013; 11: 234-44.
- 10 Soria JM et al. A quantitative-trait locus in the human factor XII gene influences both plasma factor XII levels and susceptibility to thrombotic disease. *Am J Hum Genet.* 2002; 70(3):567-574.

Description
pathologie

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Exemple

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
FAQ

Bibliographie

Bibliographie

- 11 Soria JM et al. Economic evaluation of Thromboincode, a genetic platform for the assessment of venous thromboembolism (VTE) risk, in patients with a pattern of VTE or a condition that suggests a heredity component. ISPOR 15th Annual European Congress, Berlin 3th-7th November 2012 .
- 12 Souto JC et al. Genetic susceptibility to thrombosis and its relationship to physiological risk factor: The GAIT study. Am J. Hum. Genet 2000; 67:1452-1459.
- 13 Vila J et al. Venous thromboembolism risk assessment with a multilocus genetic risk score. XX Congress of ITSH. Amsterdam 29th June-4th July 2013. PA 4.18.1
- 14 Vila J et al. Family history and genetic risk scores in venous thromboembolism risk assessment. 46th Nord Coagulation Congress. Tromso 6th-7th. Poster
- 15 Wells PS et al. Factor XIII Val34Leu Variant is protective against venous thromboembolism: A huge review and meta-analysis. Am J Epidemiol 2006; 164: 101-109.

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CRYOPEP

83 rue Yves Montand

34080 Montpellier

Tel: +33(0)4 67 10 71 20

Fax: +33(0)4 67 10 71 21

www.cryopep.fr