

TECHNOVIEW UFH

REF	5090070	TECHNOVIEW UFH CAL Set	5 x 1 mL
REF	5090071	TECHNOVIEW UFH CAL Kit	2 x 1 mL
REF	5090072	TECHNOVIEW UFH CONT L	5 x 1 mL
REF	5090074	TECHNOVIEW UFH CONT H	5 x 1 mL

symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Klíčova slova / Značenje simbola

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PRODUCT DESCRIPTION
INTENDED USE

TECHNOVIEW UFH Calibrator Set is a set of calibration plasmas and TECHNOVIEW UFH Control Plasma is a set of low and high control plasmas for calibration and the quality control of unfractionated Heparin (UFH) measurements titrated and optimised using anti-Xa assays.

COMPOSITION

TECHNOVIEW UFH Calibrator and Control plasma are prepared from citrated plasmas supplemented with different concentrations of unfractionated Heparin. The plasma contains stabilizers but no bactericide additives.

TECHNOVIEW UFH Calibrator Set 5x1 mL [REF] 5090070 contains:

mL	Reagent
1	Calibrator 1 human plasma, freeze dried, without any additive of UFH (CAL 1 at 0 IU/mL)
1	Calibrator 2 human plasma, freeze dried, supplemented with UFH (CAL 2 at ~ 0.35 IU/mL)
1	Calibrator 3 human plasma, freeze dried, supplemented with UFH (CAL 3 at ~ 0.70 IU/mL)
1	Calibrator 4 human plasma, freeze dried, supplemented with UFH (CAL 4 at ~ 1.05 IU/mL)
1	Calibrator 5 human plasma, freeze dried, supplemented with UFH (CAL 5 at ~ 1.40 IU/mL)

TECHNOVIEW UFH Calibrator Kit 2x1 mL [REF] 5090071 contains:

mL	Reagent
1	Calibrator 1 human plasma, freeze dried, without any additive of UFH (CAL 1 at 0 IU/mL)
1	Calibrator 5 human plasma, freeze dried, supplemented with UFH (CAL 5 at ~ 1.40 IU/mL)

TECHNOVIEW UFH Control Low 5x1 mL [REF] 5090072 containing 5 vials:

mL	Reagent
1	Control Low human plasma, freeze dried, supplemented with UFH (CONT L at ~ 0.20 IU/mL)

TECHNOVIEW UFH Control High 5x1 mL [REF] 5090074 containing 5 vials:

mL	Reagent
1	Control High human plasma, freeze dried, supplemented with UFH (CONT H at ~ 0.50 IU/mL)

MATERIAL REQUIRED (not supplied with the kit)

- Pipettes
- Distilled water

WARNING AND PRECAUTIONS

- IVD for *in vitro* diagnostic use
- All blood and plasma samples and products have to be regarded as potentially infectious and handled with appropriate care and in compliance with the biosafety regulations in force and must be disposed of in the same way as hospital waste.
- Each single donor plasma and each lot of Control and Calibration plasma are tested and found negative for Hb_sAg, HIV 1/2 Ab and HCV Ab. However, universal precautions (treating all human source materials as if potentially infectious) should be exercised.

STABILITY AND STORAGE

The expiry date printed on the labels applies to storage of the unopened bottles at +2...8 °C. Stability after reconstitution:

RT*	+2...8 °C	-20 °C
48 hours	7 days	1 month

The vials can be only frozen once.

Upon storage, caps should be screwed tightly.

*=room temperature ¹

TEST PROCEDURE

PREPARATION AND PERFORMANCE OF THE TEST

- Open the vial carefully and reconstitute the contents in 1 mL of distilled water by carefully rotating the vial until the product is completely reconstituted (avoid frothing).
- Allow the reconstituted TECHNOVIEW UFH Heparin Calibrator and Control plasma to stand for 10 minutes at room temperature before use ¹. Invert to mix before use (avoid foaming).
- Treat the reconstituted plasma as a citrated patient sample according to the instructions of the respective test.

APPLICATIONS FOR INSTRUMENTS

Application sheets are available from Technoclone or your local distributor upon request.

ANALYSES RESULTS

EVALUATION OF THE RESULTS

The figures in the batch table are only applicable to the indicated lot number of TECHNOVIEW UFH Calibrator and Control plasma and the method has to adhere strictly as prescribed for the respective reagent.

STANDARDIZATION

TECHNOVIEW UFH Calibrator and Control plasma are calibrated against the 5th International Standard preparation for UFH 97/578.

LIMITATION OF THE TEST

The figures in the batch table are only applicable to the indicated lot number of TECHNOVIEW UFH Calibrator and Control plasma and the method has to adhere strictly as prescribed for the respective reagent.

The guidelines for evaluation given in this leaflet have been established assuming that a constant sensitivity of different batches is guaranteed by the manufactures of each reagent.

LITERATURE

1. Richardson J.H. and Backley W.E. Eds. Bios. Microbiological and Biomedical laboratories. US. Dept. and Human Services, Public Health Service, HHS publication (CDC) 84-8395, Washington, D.C. (1984).
2. National Committee for Clinical Laboratory St Specifications for reagent water used in the clinical lab NCCLS Approved Standard: ASC-3.
3. Westgard J.O., Barry P.L. Cost effective Quality Managing for managing the quality and the Productivity of Analytical AACC Press (1986).

¹ For standardization a reconstitution time of 30 minutes is recommended.

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