# TECHNOZYM<sup>®</sup> PAI-1 Actibind<sup>®</sup> ELISA

REF	TC16075	TECHNOZYM <sup>®</sup> PAI-1 Actibind <sup>®</sup> ELISA	
REF	TC16077	TECHNOZYM <sup>®</sup> PAI-1 Actibind <sup>®</sup> Calibrator Set	5 x 0.2 mL
REF	TC16079	TECHNOZYM <sup>®</sup> PAI-1 Actibind <sup>®</sup> Control Set	2 x 0.2 mL

## symbols key / Symbolschlüssel / interpretazione dei simboli / explicación de símbolos / explicação dos símbolos / clé des symboles / Symbolnyckel / symbolforklaring / Tegnforklaring / Κλειδί συμβόλων / Използвани символи / символы / Ključova slova / Značenje simbola

	manufacturer / Hersteller / fabbricante / fabricante / fabricante / fabricant / Tillverkaren / Fabrikanten / Produzent / Κατασκευαστής / Производитель / Производител / Производител / výrobce / Proizvođač		expiry date / Verfallsdatum / data di scadenza / fecha de caducidad / data de validade / date d'expiration / utgångsdatum / udløbsdato / Utløpsdato / Ημερομηνία λήξης / срок на годност / datum expirace/ срок годности / datum expirace / Rok trajanja
	storage temperature / Lagertemperatur / temperatura di conservazione / temperatura de conservación / temperatura de conservação / température de stockage / lagringstemperatur / opbevaringstemperatur / Oppbevaringstemperatur / θερμοκρασία αποθήκευσης / съхранение на / teplota skladování / температура хранения / teplota skladování / Temperatura lagerovanja		consult instructions for use / Gebrauchsanweisung beachten / consultare le istruzioni per l'uso / consulte las instrucciones de uso / consultar o manual de instruções / instruction d'utilisation / se användarinstruktioner / følg brugsvejledning / Følg bruksanvisningen / συμβουλευθείτε τις οδηγίες για τη χρήση / прочетете инструкцията за работа / potfeba fidiť se instrukcemi / перед использованием читайте инструкцию / sledujte návod k použití / Pročitaj upustvo pre upotrebe
	CE-mark / CE-Kennzeichnung / marchio CE / marca de CE / Simbolo da CE / marquage CE / CE-märkning / CE-mærket / CE-merke / CE-σημάδι / CE марка / CE-označení / маркировка CE / značka CE / CE-marka		determinations / Bestimmungen / determinazioni / determinaciones / determinações / déterminations / bestämmingar / bestemmelser / Bestemmelser / προσδιορισμοί / брой тестове / stanovení / определний / роčet stanovení / Definicija
<b>AQUA</b>	distilled water / destilliertes Wasser / acqua distillata / água destilada / água destilada / eau distillée / destillerat vatten / destilleret vand / Destillert vann / απεσταγμένο νερό / децтилирана вода / destilovaná voda / дистиллированная вода / destilovaná voda / Serija	<b>LOT</b>	lot / Charge / lotto / lote / lote / lot / sats / serie / Parti / парти / партида номер / šarže / lot / šarže / in vitro diagnostika
<b>BUF</b>	Reaction buffer / Reaktionspuffer / tampone di reazione / tampón de reacció / Tampão de reação / tampon de réaction / Reaktionsbuffert / Reaktionsbuffer / Reaktionsbuffer / διάλυμα αντίδρασης / Реакционен буфер / Рабочий буферный раствор / Reakčni pufr / Reakcioni pufer	<b>MTP</b>	microtiter plate / Mikrotiterplatte / placa microtiter / microplaca / microplaca / microplaques sensibilisées / Mikrotiterplatta / Mikrotiterplade / mikrotiterplate / πλάκα μικροτιτλοδότησης / Микротитрѐна плака / Микропланшет / Mikrotitrační destička / Mikrotitracione ploče
<b>CAL</b>	Calibrator / Kalibrator / Calibratore / calibrador / calibrador / calibreur / Kalibrator / Kalibrator / Kalibrator / Βαθμονομητής / Калибратор / калибратор / kalibrátor / Kalibrator	<b>REF</b>	catalogue number / Katalognummer / numero di catalogo / número de catálogo / número de referência / réf. de catalogue / katalognummer / Katalognummer / αριθμός καταλόγων / каталожен номер / katalogové číslo / каталожный номер / katalogové číslo / Kataloški broj
<b>CONJ</b>	Conjugate / Konjugat / Coniugato / conjugado / conjugado / conjugaté / Konjugerad / Konjugat / Konjugat / συνδεδετικό / Конюгат / Коньюгат / Konjugát / Konjugat	<b>RTU</b>	ready to use / gebrauchsfertig / pronto all'uso / listo para usar / pronto a usar / prêt à l'emploi / færdig att användas / færdig til brug / klar til bruk/ έτοιμο προς χρήση / Готов за употреба / готов к использованию / k přímému použití / Razrediti ili rastvoriti
<b>CONT</b>	Control / Kontrolle / controllo / control / control / contrôle / Kontroll / Kontroll / Kontroll / διάλυμα ελέγχου / Контрол / Контрольный образец / Kontrola / Kontrola	<b>STOP</b>	stop solution / Stopplösung / Soluzione di arresto / solución de parada / solução de paragem / solution d'arrêt / Stopplösning / Stop-opløsning / Stoppløsning / διάλυμα παύσης / Стоп раствор / Стоп-раствор / Zastavovací roztok / Stop solucija
<b>DIL</b>	dilute or dissolve in / verdünnen oder lösen in / diluire o dissolvere in / diluir o dissolver / diluir ou dissolver em / diluer ou dissoudre dans / spädd eller upplös i / fortyndes eller opløses i / Fortyndes eller opløses i / αραιωση ή διάλυση σε / растворите или разрежете с / zředit anebo rozpustit v / разбавить или растворить в / naředte nebo rozpustte v / razrediti ili rastvoriti u	<b>SUB</b>	substrate / Substrat / substrato / substrato / substrato / substrat / Substrat / Substrat / Substrat / υπόστρωμα / Субстрат / Субстрат / Substrát / Substrat
<b>INC</b>	incubation buffer / Inkubationspuffer / tampone di incubazione / tampón de incubación / tampão de incubação / tampon d'incubation / Inkubationsbuffert/ Inkubationsbuffer/ Vaskebufferkonsentrat / διάλυμα επώασης / Инкубационен буфер / Буфер для инкубации / Inkubační pufr / Inkubacioni pufer	<b>WASH</b>	washing solution concentrate / Waschlösungskonzentrat / concentrado de solución de lavado / solución de lavado concentrada / tampão de lavagem concentrado / Tampon de lavage concentré / Vattenlösningsskoncentrat / Vaskeopløsningskonsentrat / vaskeløsningskonsentrat / συμπυκνωμένο διάλυμα πλύσης / Концентриран миеш разтвор / Концентрат промывочного раствора / Koncentrát promývacieho roztoku / Koncentrat solucije za ispiranje
<b>IVD</b>	in vitro diagnostic use / in vitro Diagnostikum / diagnostico in vitro / diagnóstico en vitro / diagnóstico in vitro / diagnostic in vitro / för in vitro diagnostik / in vitro diagnostik / in vitro diagnostisk bruk / рηση διαγνωστικής εντόβ σωλήνα / за ин vitro диагностика / pro in vitro diagnostiku / использовать для диагностики in vitro / diagnosticky prostfedek in vitro / Destilísana Voda		

**PRODUCT DESCRIPTION**

**INTENDED USE**

The TECHNOZYM® PAI-1 Actibind® ELISA can be used to determine active PAI-1 levels in patients with thrombotic disorders (deep vein thrombosis, myocardial infarction stroke), malignancies or septicemia.

**COMPOSITION**

1. ELISA test strips (12); with 8 wells each, precoated with tPA immobilized via a monoclonal anti tPA coating antibody; in an aluminium bag.
2. Washing buffer concentrate (PBS; pH 7.3); containing detergent; 0.01% merthiolat; 1 bottle, 80 ml.
3. Incubation buffer (PBS; pH 7.3); contains stabiliser protein; 0.05% proclin; and dye, 1 bottle, 90 ml, ready for use.
4. Calibrators (Standards) numbered; lyophilised; 1 bottle each. **Concentrations are lot-dependent; consult label on the vial.**
5. Control plasmas "low level" and "high level" for checking purposes lyophilised; 1 bottle each. **Concentrations are lot-dependent; consult the label on the vial.**
6. Conjugate monoclonal Anti-PAI-1-POX; dyed blue; 1 bottle, 0.3 ml.
7. Chromogen TMB (tetramethylbenzidine); 1 bottle, 12 ml; ready for use.
8. Stopping solution sulphuric acid 0.45 mol/l; 1 bottle 12ml; ready for use.
9. Adhesive film: for ELISA test strips (2).

**MATERIAL REQUIRED** (but not supplied with the kit)

1. Distilled water
2. Test tubes for diluting standard and samples
3. Measuring cylinder (1000 ml)
4. Precision pipettes (10, 100 and 1000 µl)
5. Variable pipette (1000 µl)
6. Multichannel and/or dispensing pipettes (100 and 200 µl)
7. ELISA washer or multichannel pipette
8. ELISA reader with 450 nm filter, with a 620 nm reference filter if available.
9. Incubator (37 °C)

**WARNING AND PRECAUTIONS**

- All human blood or plasma products as well as samples must be considered as potentially infectious. They have to be handled with appropriate care and in strict observance of safety regulations. The rules pertaining to disposal are the same as applied to disposing hospital waste.
- Calibrators and control plasmas are made from human blood and any individual plasma involved in the procedure is HBsAg, HIV 1/2 Ab and HCV-Ab-negative (see labels on kit and/or bottles).
- Stopping solution (sulphuric acid) may irritate the skin. Should acid get into your eyes, wash out immediately with water and consult a doctor.
- The reagents sometimes contain preserving agents (merthiolate). Beware of swallowing! Avoid contact with skin or mucous membranes.

**STABILITY AND STORAGE**

The expiry date printed on the labels applies to storage of the unopened bottles at + 2...8 °C.

Stability after reconstitution/opening:

Material/ Reagent	State	Storage	Stability
Calibrators, controls plasmas	after reconstitution	-20 °C	6 months
ELISA test strip	after opening	2 ... 8 °C with adhesive film in plastic bag with drying agent	expiry date
Washing buffer concentrate	after opening	2 ...8 °C	6 months
Washing buffer	1+11.5 dilution of concentrate	2 ... 8 °C	3 weeks
Incubation buffer	after opening	2 ... 8 °C	2 months
Conjugate	after opening	2 ... 8 °C	6 months
	working solution	room temperature	60 minutes
Chromogen TMB	after opening	2 ... 8 °C	expiry date

**TEST PROCEDURE**

**PREPARATION OF SAMPLES**

Sample material: plasma  
 It is highly recommended to use commercially available collecting tubes which contain platelet stabilizing agents e.g. CTAD (Greiner). 90% of PAI-1 antigen is contained in the platelets so it is essential to ensure during sample collection that the platelets are not damaged which would result in elevated plasma levels. Citrated or EDTA plasmas can be used. Centrifuge for 15 minutes at a minimum of 2500g (DIN 58905). The plasma sample may be stored for 3 hours at room temperature; otherwise the sample ought to be frozen immediately after centrifugation at -20°C or below. Plasmas are stable for 6 months at -20°C.  
 Serum samples should not be used as they show high PAI-1 levels which is related to the platelet PAI-1 content. Cell supernatants and tumour extracts can be used, but this ELISA test has been optimized for plasma samples, therefore other dilution factors would have to be used accordingly.

**PREPARATION OF REAGENT**

1. Before starting the test, all the required components are to be brought to room temperature.
2. Preparing the washing buffer: Dilute 1 part by volume washing buffer concentrate with 11.5 parts by volume distilled water (1+11.5). Mix well! (Diluted washing buffer concentrate = washing buffer). There may be crystalline precipitations which will dissolve at 37°C within 10 minutes.
3. Reconstituting calibrators and control plasmas: Calibrators and control plasmas are reconstituted with 200 µl distilled water and mixed for 10 seconds after a reconstitution time of 15 minutes (vortex mixer). Reconstituted components are clear to slightly turbid.
4. Preparing the conjugate working solution (1+50): Dilute 1 part by volume conjugate with 50 parts by volume incubation buffer.

**For 8 test wells: Mix 20 µl conjugate with 1000 µl incubation buffer**

**PERFORMANCE OF THE TEST**

<b>WASHING</b> (reference 1,3,4)	washing buffer	5 x 200 µl
<b>SAMPLE AND CONJUGATE REACTION</b> (reference 1,2)	Calibrators, control plasmas and samples	25 µl
	conjugate working solution	75µl
	Pipette in wells, cover test, strips with film, incubate at 37°C	45 minutes
<b>WASHING</b> (reference 1,3,4)	washing buffer	5 x 200 µl
<b>SUBSTRATE REACTION</b> (reference 1,2)	pipette substrate solution into test wells	100 µl
	Cover test strips with film, incubate at RT (20 – 25°C)	15 minutes
<b>STOPPING</b> (reference 1,2)	pipette stopping solution into wells	100 µl
<b>MEASURING</b> (reference 5)	ELISA-Reader, 450 nm	shake 10 sec., measure within 10 minutes

**References**

1. Reagents of different lots must not be combined
2. Precision and performance, among others, essentially depend on the following factors:
  - Thorough mixing of all substances used for dilution
  - Test calibrators, controls and samples in duplicates.
  - Incubation to be done at correct temperatures
  - Strict observance of the order of pipetting and of the time element as indicated:
  - The time for sample incubation, conjugate and substrate reaction as indicated starts after pipetting the last sample. Incubation times should not vary by more than ±10%.
  - During sample incubation and conjugate reaction, the time for pipetting the diluted calibrators/samples/control plasmas and/or conjugate solutions must not exceed 60 seconds per ELISA test strip (8 wells).
  - During substrate reaction and at stopping, the time needed for pipetting the substrate and/or the stopping solution must not exceed 10 seconds per ELISA test strip. Short pipetting times may be secured by using multichannel- and dispensing pipettes.
3. Label/number strips with a water resistant pen in case the strips accidentally fall out of the frame during testing.
4. After the last washing, wells must be aspirated thoroughly, turned upside down and positioned on a blotting paper; by gentle tapping, the last remnants must be removed.
5. Measuring the difference in wave lengths at 450 and 620 nm or at 450 and 690 nm, the precision of the test is increased.

**LIMITATION OF THE TEST**

Samples which fall higher than the top calibrator standard must be retested at a higher dilution as a hook dose response may occur above.

**ANALYSIS RESULTS**


**CALCULATION OF THE RESULTS**

Setting up a reference curve: X axis: Concentration active PAI-1 U/ml  
 Y axis: Extinction  
 Graph plot is linear-linear with a linear or point to point fit

Assessment of reference curve

- The extinction coefficient of the highest calibrator should be between 1.0 and 2.5.
- The validity of the test may be checked on the basis of the calculated control values.

Example of standard curve.

Measuring concentration of samples

- Read off the concentration from the reference curve.
- If there are samples with extinction coefficients higher than that of the highest point on the curve, they have to be prediluted with incubation buffer (1:2 or 1:4). The measured concentration then has to be multiplied with the dilution factor 2 or 4.

**REFERENCE RANGE**

Normal Plasma levels range from 1-7U/ml. It is recommended that individual laboratories establish their own normal range. Active PAI-1 levels above 20U/ml may indicate reduced fibrinolytic capacity and, thus, increased thrombotic tendency. Measures should be taken to reduce the risk of thrombosis in individuals with elevated PAI-1 plasma levels. The TECHNOZYM® PAI-1 Actibind® ELISA exclusively measures free, active PAI-1 and is not affected by other forms of PAI-1 or other plasminogen activator inhibitors.

**STANDARDISATION**

The calibration material used is the WHO International Standard for Plasminogen Activator Inhibitor (PAI-1).

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