

CERTIFICAT D'ANALYSES

pNAPEP 1025
Lyophilisat

Référence
61031025

Lot
180702

Péremption
2024/04

Usage RECHERCHE Exclusivement

RUO

2°C / 8°C

Conditions de stockage : Entre +2/+8°C à l'abri de la lumière et de l'humidité
Type : Substrat du facteur Xa
Séquence : CH₃SO₂-(D)Leu-Gly-Arg-pNA, AcOH
Poids moléculaire : 602,6 g/mol

	<i>Spécification</i>	Résultat	Conformité
Aspect	<i>Poudre blanche à jaune pâle</i>	Poudre blanche	Conforme
pH	<i>5 < < 8</i>	5,07	Conforme
pNA libre	<i>< 0,5 %</i>	0,34 %	Conforme
Teneur en peptide	<i>≥ 25 mg</i>	26,66 mg	Conforme
Reconstitution	<i>2 mmol/L</i>	22,12 mL d'eau distillée	

Contrôle Qualité



CERTIFICATE OF ANALYSIS

pNAPEP 1025
Lyophilisate

Catalog number
61031025

Lot
180702

Expiry
2024/04

RESEARCH USE ONLY

RUO

2°C / 8°C

Storage conditions : +2/+8°C avoid light exposure and moisture
Type : Activated factor X substrate
Identity : CH3SO2-(D)Leu-Gly-Arg-pNA, AcOH
Molecular weight : 602,6 g/mol

	<u>Specification</u>	Result	Conformity
Appearance	<i>White to pale yellow powder</i>	White powder	Conform
pH	<i>5 < < 8</i>	5,07	Conform
Free pNA	<i>< 0,5 %</i>	0,34 %	Conform
Quantity /vial	<i>≥ 25 mg</i>	26,66 mg	Conform
Reconstitution	<i>2 mmol/L</i>	22,12 mL of distilled water	

Quality Control

